

OBJECTIVE: Preparation of competent cells from *E. coli* DH5 α cells

Chemicals required: MgCl₂, CaCl₂, glycerol, luria agar (HiMedia), luria broth (HiMedia), ddH₂O/MilliQ.

Plastic ware/ glassware: measuring cylinder, autoclaved micro centrifuge tubes, 1.0 ml tips, 0.02ml-0.2 ml tips, enzyme tips, 100ml/250 ml flasks, reagent bottles, Test tubes, OakRidge tubes, Uvettes.

Equipment required: refrigerated centrifuge, vortex, UV spectrophotometer,

Recipe for solutions:

Luria broth: 2.0 g in 100 ml (autoclaved)

Solution I (100 ml)

S.No.	Chemical(s)	Concentration	Weight
1.	MgCl ₂	80 mM	1.642 g
2	CaCl ₂	20 mM	0.294 g

Solution II (100 ml)

S.No.	Chemical(s)	Concentration	Weight
1.	CaCl ₂	0.1 M	1.4702 g

Solution III (100 ml)

S.No.	Chemical(s)	Concentration	Weight
1.	CaCl ₂	0.1 M	1.4702 g
2	Glycerol	15%	

Solutions were prepared and autoclaved

PROTOCOL:

Inoculate a single colony of *E. coli* DH5 α in 5ml LB (primary inoculum) and incubate at 37°C, 180 rpm overnight.

Transfer 1 percent of primary inoculum in 100 ml LB and grow at 37°C, 180 rpm until O.D_{600nm} reaches ~0.5-0.6.

After attaining appropriate absorbance keep the culture on ice for 20 min.

Transfer the culture in pre-chilled OakRidge tube and harvest cells at 4000 rpm for 10 min, 4°C.

Discard supernatant and resuspend the pellet in 30 ml of chilled Solution I

Recover cells by centrifugation at 4,000 rpm for 10 min at 4°C.

Discard supernatant and resuspend the pellet in 2-4ml of solution II.

Recover cells by centrifugation at 4,000 rpm, for 10 min at 4°C.

Discard supernatant without disturbing pellet, and resuspend in 0.8 ml-1.0 ml of solution III.

Prepare aliquots of 100µl and store at -80°C.